

The most compact intelligent ICCD picos-fast camera system (digital or analog output)



Features

- Fastest shutter down to 0.2ns
- Gating times from 0.2ns \rightarrow DC
- 18mm standard, 25 mm image intensifier optional
- Spectral Sensitivity of photocathode from UV NIR (110 1300nm)
- Multiple Exposures
- Single Photon detection
- Analog (EIA/CCIR) or digital output (12, or 14 Bit with Standard or High Resolution)
- Customized f/0.8 distortion free lens coupling between image intensifier and CCD
- High Dynamic Range: 14 Bit (21 Bit with 4 Spec E spectroscopy software)
- Effortless Image/Data storage and retrieval
- Single Trigger Discriminator integrated
- Internal or external trigger
- Free Terminal software and printed manual

SPECIFICATIONS

Unique Features	4 Picos
Shortest gating time	Shortest gating time 200ps
Intensifier output coupling	customized distortion free f/0.8 relay lens
Optical input	c-mount (standard), Nikon F-mount (optional)
Optical input window	fused silica
Sensitivity corresponding to conventional film	1×10^{12} fc = 50 x 5 ⁷ ASA
Multiple exposures, "dead" time between exposures	any sequence 0.3µm
Gate repetition rate	3.3MHz burst, 200kHz continuous
TTL pulse (incorporated)	standard
Remote control	RS 232 (digital set up)
CE certified	ves

Adaptation to in situ light level variations by internal digital programmable brightness control of the electronic shutter action, with exact reproducible digital setting of delay and exposure times.

Very high system integration permits small physical size of the total unit (all in one) — even very difficult surveillance jobs are mastered easily by remote control. Programmable control parameter entry via RS 232 (digital set up), remote control software included.

Distortion free imaging due to advanced proximity focused MCP (Micro Channel Plate) image intensifier and use of highest quality CCD array for best sensitivity and resolution.

14 Bit High Dynamic Range (theoretical limit 16 bit), with 4 Spec E Spectroscopy PC Software up to 21 Bit/ Spectrum with all lines integrated.

Automatic Exposure Control and Automatic Gain, for unattended operation under greatly varying light conditions (optional). Effortless Image/Data storage and retrieve via system interface RS 170 to an external VCR (standard), frame grabber or optical disk (optional).

Switches and Connections

No.	Item	Function
А	Switch	Power Switch ON/OFF
В	Socket	Power Supply Socket (12V)
1	Video	Camera output signal (RS170
2	Busy	Synchronization Signal (TTL)
3	V _{Init}	Asynchronous reset of CCD of
4	F_{Sync}	TTL output for synchronization
5	-Trig	Trigger input, negative edge 1
6	+Trig	Trigger input, positive edge T
7	IntGtP	Output of internal time delay generator
8	ExtGtP	Input for control of HV MCP p



For control by internal time/delay generator, 7 and 8 are shorted internally. Camera will free run when 4 to 5, and 7 to 8 are internally connected. These settings are both under RS 232 control. Camera may be externally driven through 8 by external pulse delay generator when 7 and 8 are disconnected. Pulse monitoring is provided by 7 and camera master sync ouput is available at 4.

Rear view of 4 Picos with legend



ICCD - pSec



www.toptical.com.tw 02-2346-1510

- or CCIR)
- e.g. Frame Grabber
- amera
- n purposes
- TTL
- TL
- gate pulse
- pulse (TTL)

oscopes	
 Imaging 	CCD Cameras
Communications	Semiconductors ·
Lighting	Solar Cells ·
Instruments	Tests ·
Sensors	Detection •
Mechanics · Positioning	Components ·
Light Sources	Lasers ·



Image Intensifier

Image Intensifier specifications

Image intensifier type (proximity focused MCP)	single stage (standard), dual stage (optional)
Phosphor material	P43 (standard), P46 (optinal)
Image intensifier diameter (mm)	18mm standard, 25mm optional
Image area of the relay lens	25mm MCP: 20 x 15mm, 18mm MCP: 14.4 x 10.8mm
Wavelength range, subject to window design	165 – 920nm (quartz, standard), 110 – 1300nm (optional)
Spectral Sensitivity of MCP (nm)	1101300nm, depends on the type of the photocathode
Quantum Efficiency (Q.E.) (see curves below)	depends on the type of the MCP, up to 35%
Gain (4k steps) (V _{MCP} =1000V) control via RS 232 digital setup	single stage MCP: 4 x 10 ⁴ dual stage MCP: 4 x 10 ⁶
Coupling phosphor (MCP \rightarrow CCD)	customized 6 element f/0.8 relay lens No distortion! No vignetting! No pin cushion!

Spectral Sensitivity of Photocathodes (Wavelength in nm)

Standard 18 mm		Optional 25 mm			
S20UV	В	approx. 165 - 820nm	Solar Blind (CsTe)	G	approx. 180 - 340nm
S25 IR (Super S25)	н	approx. 350 - 920nm	Bialkali		approx. 165 - 600nm
Optional 18 mm		Enhanced S20	D	approx. 165 - 820nm	
S20 UV (MgF2)	А	approx. 110 - 820nm	Enhanced S25 (glass)	I	approx. 270 - 900nm
Broadband	J	approx. 190 - 920nm	Wideband S25 WB	к	approx. 200 - 900nm
Standard 25 mm		S1	Е	approx. 700 - 1300nm	
S20	С	approx. 165 - 820nm	Enhanced S20	D	approx. 165 - 820nm
S25	F	approx. 200 - 840nm			



Deviations of up to $\pm 25\%$ from the above typical spectral sensitivity curves are possible. The position of the curves can vary ± 20 nm. The input window material limits the spectral response of the photocathode in the shorter wavelengths. The window materials and their transmission limits are: quartz (165nm), MgF2 (110nm).

Image intensifier and shutter (schematic)



CCD-Video Unit

CCD Video Chin	Analog Out	tput	Standard	High
CCD video Chip	USA, Japan	Elsewhere	Resolution CCD	Resolution CCD
Analog or digital output	analog EIA (RS 170)	analog CCIR	digital 12 or 14 Bit	digital 12 or 14 bit
Resolution (pixel)	768 x 494	752 x 582	780 x 580	1360 x 1024 (1x1) 580 x 512 (2x2) 456 x 342 (ROI)
Pixel size (µm)	8.4 x 9.8	8.6 x 8.3	8.3 x 8.3	4.7 x 4.7
Imaging frequency (analog) Frame rate (digital)	30/60Hz	25/50Hz	12bit/14bit: 33.8/60.8/67.0 fps	12bit/14bit: 10.6/17.9/20.9 fps
Dynamic Range A/D (Bit)	14 Bit, up to 21 Bit (wit	h 4 Spec E spec	troscopy software)	
Video Gain	025dB, automatic adjustable through computer RS 232 inter	or manually face	Full Frame, ROI: 02 Binning: 025dB	0dB
Binning vertical (pixel)	Software		1 (full frame) ,2 (Binn	ing) pixel, ROI
Binning horizontal (pixel)	Software		1 (full frame) ,2 (Binn	ing) pixel, ROI
Image Sensor	ICXAL			
Chip Readout	Correlated double sam	pling, dark curre	nt corrected	
Output	$1V_{PP}$ (75 Ω), composite	e video, RS 170/I	EIA, CCIR or VGA	
Scan Mode	Field/Frame, selectable through computer RS 232 interface. ICCD camera 4 Picos or 4 Quik E may be Genlocked or supply Fsync-pulse to operate as master clock.			
Gamma	1 or .45, selectable three	ough computer F	RS 232 interface	
Internal Synchronization	Free run mode			
External Synchronization	by negative edge TTL	pulse (Vinit)		
Cooling of CCD optional	Regulated cooling of C of 10 for exposure ti condensation; eliminat	CCD camera uni mes above 100 es need for vacu	t to 14°C to minimize on ms. Provides single um or special nitrogen	dark current by a factor photon sensitivity. No atmosphere.

t

ICCD - pSec



www.toptical.com.tw 02-2346-1510 t



4 Picos control window

			5		Set Mode
Unit		[us] 💌			Single
Sequence	0,50,200,50,200,1	00			
Varning Tim	e exceeds 100 e	R Ac	cept ?		
MCP Gain	Voltage 720	_			
Scan Mode	Video A	cq.	Frame Gr	abber	
Frame	G Stande G Stande G Stande Stande	brd	@ None		Bend h
C Field	C CCD Ir	it.	C On Bos	urd	
SyncOutput	Trig. Source	Trigger	Mode	Gate Ctrl	Connect
C Vertical	@ Fsync	CDirect	Trigger	@ IntGt	Initialize
e Odd	C -Trig	@ STD	1	CExtGt	
Gamma	Video Gain		Start Optio	n	Port
0.45	C Auto		C Cold		COMI
61	G Eined J	C	Cilline	G Auto	

Shutter control

The advanced, digitally controlled shutter delay feature is the perfect match for your laser, range gating, flow analysis, or many other high speed applications. It is operational in 'scattered light' environments, underwater or for highest speed multi-instrument sequential image acquisition. Multiple direct images with a repetition/ delay time setting as short as 0.3µs can be synchronized with ease to any external TTL source.

Internal exposure control	4 Picos
Time (t_s) and delay (t_D) of the gate pulse, or multiple exposure with CPU internally digital programmable	t_s = 200ps 80s, min. steps 10ps t_D = 0 80s, min. increments 10ps
Trigger propagation delay	<65ns, less than 10ps jitter
Initializing	-Trig, +Trig, or FSync
Multiple Exposure	Any sequence, 0.3µs "dead time" between exposures

External exposure control

4 Picos

Control of the camera internal Pulse E amplifier via ExtGtP (TTL Pulse) input: Shutter continuous from:

<45ns, no jitter

 $t_{\rm S}$ = 200ps ... , $t_{\rm D}$,

t, t_D determined by external device

Trigger propagation delay

External exposure control

Pixel by pixel exposure analysis providing automatic light level control by instantaneous adjustment of camera shutter speed and intensifier gain for very wide range of lighting conditions (up to 12 orders of magnitude). 200ps ... 15 (18) ms, Shutter time and MCP gain automatically adjusted, in response to scene illumination.



Analog shutter control (schematic)

Mechanical & Environmental Data

Mechanical & Environmental Data, Power Requirements

Camera dimensions, without lens (mm, inch)	248 9 ³ /4
Camera weight (all in one) (kg / lb)	3kg
Camera mount (at the bottom plate of the camera)	1⁄4"
Operating Humidity (%)	25.
Operating temperature (°C / °F)	0°C
Performance specification	10°
Operating limits	-10
Shock and Vibration	60g
Voltage	12



Ŕ



www.toptical.com.tw 02-2346-1510 toptical@ms2



 $8 \times 110 \times 135$ mm (l x w x h) $4 \times 4^{3}/8 \times 5^{4}/16^{\circ}$ (l x w x h)

g / 6.6lb

x 20 and M8 mounting hole

.95%, non condensing

C – 50°C / 32°F – 122°F

°C – 40°C / 50°F – 104°F

°C – 50°C / 14°F – 122°F

g accel. shock, 7g Vibration (11 – 200Hz)

V +5%/-2%

Spectroscopes	
 Imaging 	CCD Cameras
Communications	Semiconductors ·
Lighting	Solar Cells ·
Instruments	Tests ·
Sensors	Detection .
Mechanics · Positioning	Components ·
Light Sources	Lasers ·
2 -	89



4 Picos/... intelligent ICCD Camera Family

analog or digital output



5ns 6ns 7ns





10ns

Analysis of plasma development with 4 Picos camera: single beam set up. Plasma images of water droplets . Shutter time: 200ps. Laser Pulse: 50mJ, 8ns

© 2003 with courtesy of Ch. Janzen, Fraunhofer Institut Aachen, Germany





Multiple exposures (Mouse trap)

Multiple 16 bit exposures of Xe flash tube, single discharge, 10 exps during 120µs, continuously swept by scanning mirror.

Ŕ





Applications

Physical Sciences

Plasma temperature and density analysis Plasma flow analysis Combustion analysis Synchrotron radiation Laser induced fluorescence

Engineering Research

Particle Tracking Velocimetry (PTV) Particle Imaging Velocimetry (PIV) Automative Fuel Injection Spray analysis Wind tunnel studies Stress analysis of ceramics materials

Biological Sciences

Cancer research Fundus imaging spectroscopy X-ray detection Luminescence Time resolved fluorescence

High Speed Imaging

Dynamic Schlieren Phenomena Shock tubes Range gating

Low Light Imaging

Thomson Scattering Raman Spectroscopy Glow Discharge Spectroscopy Semiconductor failure analysis





Nitrogen-laser-based system for oral cancer diagonis developed at Center for Advanced Technology, Indore, India © with courtesy of CAT, Indore, India

4 Quik E

The most compact intelligent ICCD nano-fast camera system (digital or analog output)

Features

- Fastest shutter down to 1.2ns typically
- Gating times from 1.2ns typically \rightarrow DC
- 18 or 25 mm Image intensifier (to be specified)
- Spectral Sensitivity of photocathode from UV NIR (110 1300nm)
- Multiple Exposures
- Single Photon detection
- Analog (EIA/CCIR) or digital output (12, or 14 Bit Standard or High Resolution)
- Customized f/0.8 distortion free lens coupling between image intensifier and CCD
- High Dynamic Range: 14 Bit (21 Bit with 4 Spec E spectroscopy software)
- · Effortless Image/Data storage and retrieval
- Single Trigger Discriminator integrated
- Internal or external trigger
- · Free Terminal software and printed manual

SPECIFICATIONS

Unique Features	4 Picos
Shortest gating time	Shortest gati
Intensifier output coupling	customized of
Optical input	c-mount (sta
Optical input window	fused silica
Sensitivity corresponding to conventional film	1 x 10 ¹² fc = 5
Multiple exposures, "dead" time between exposures	any sequenc
Gate repetition rate	3.3MHz burs
TTL pulse (incorporated)	standard
Remote control	RS 232 (digit
CE certified	yes

*Actual performance may vary with manufacturer and specimen, please ask for details.

Ż

ICCD

nSec



0

www.toptical.com.tw 02-2346-1510



ing time 200ps

distortion free f/0.8 relay lens

ndard), Nikon F-mount (optional)

50 x 5⁷ ASA

e 0.3µm

st, 200kHz continuous

tal set up)

pectroscopes	
 Imaging 	CCD Cameras
Communications	Semiconductors ·
Lighting	Solar Cells ·
Instruments	Tests ·
Sensors	Detection .
Mechanics · Positioning	Components ·
Light Sources	Lasers .
2 -	93

Adaptation to in situ light level variations by internal digital programmable brightness control of the electronic shutter action, with exact reproducible digital setting of delay and exposure times.

Very high system integration permits small physical size of the total unit (all in one) — even very difficult surveillance jobs are mastered easily by remote control.

Programmable control parameter entry via RS 232 (digital set up), remote control software included.

Distortion free imaging due to advanced proximity focused MCP (Micro Channel Plate) image intensifier and use of highest quality CCD array for best sensitivity and resolution.

14 Bit High Dynamic Range (theoretical limit 16 bit), with 4 Spec E Spectroscopy PC Software up to 21 Bit/ Spectrum with all lines integrated.

Automatic Exposure Control and Automatic Gain, for unattended operation under greatly varying light conditions (optional). Effortless Image/Data storage and retrieve via system interface RS 170 to an external VCR (standard), frame grabber or optical disk (optional).

Switches and Connections

No.	ltem	Function
А	Switch	Power Switch ON/OFF
В	Socket	Power Supply Socket (12V)
1	Video	Camera output signal (RS170 or CCIR)
2	Busy	Synchronization Signal (TTL), e.g. Frame Grabber
3	V _{Init}	Asynchronous reset of CCD camera
4	F_{Sync}	TTL output for synchronization purposes
5	-Trig	Trigger input, negative edge TTL
6	+Trig	Trigger input, positive edge TTL
7	IntGtP	Output of internal time delay gate pulse generator
8	ExtGtP	Input for control of HV MCP pulse (TTL)



For control by internal time/delay generator, 7 and 8 are shorted internally. Camera will free run when 4 to 5, and 7 to 8 are internally connected. These settings are both under RS 232 control. Camera may be externally driven through 8 by external pulse delay generator when 7 and 8 are disconnected. Pulse monitoring is provided by 7 and camera master sync ouput is available at 4.

Rear view of 4 Picos with legend

Image Intensifier

Image Intensifier specifications

Image intensifier type (proximity focused MCP)	single stage (
Phosphor material	P43 (standa
Image intensifier diameter (mm)	18mm stan
Image area of the relay lens	25mm MCF
Wavelength range, subject to window design	165 – 920n 110 – 1300
Spectral Sensitivity of MCP (nm)	1101300n depends or
Quantum Efficiency (Q.E.) (see curves below)	depends or
Gain (4k steps) (V _{MCP} =1000V) control via RS 232 digital setup	single stage dual stage l
Coupling phosphor (MCP \rightarrow CCD)	customized No distortio

Spectral Sensitivity of Photocathodes (Wavelength in nm)

Standard 18 mm			Optional 25 mm		
S20UV	В	approx. 165 - 820nm	Solar Blind (CsTe)	G	approx. 180 - 340nm
S25 IR (Super S25)	Н	approx. 350 - 920nm	Bialkali		approx. 165 - 600nm
Ор	tiona	ıl 18 mm	Enhanced S20	D	approx. 165 - 820nm
S20 UV (MgF2)	А	approx. 110 - 820nm	Enhanced S25 (glass)	I	approx. 270 - 900nm
Broadband	J	approx. 190 - 920nm	Wideband S25 WB	К	approx. 200 - 900nm
Standard 25 mm		S1	Е	approx. 700 - 1300nm	
S20	С	approx. 165 - 820nm	Enhanced S20	D	approx. 165 - 820nm
S25	F	approx. 200 - 840nm			

Deviations of up to ±25% from the above typical spectral sensitivity curves are possible. The position of the curves can vary ±20nm. The input window material limits the spectral response of the photocathode in the shorter wavelengths. The window materials and their transmission limits are: quartz (165nm), MgF2 (110nm).



CCD Cameras Imaging

Communications Semiconductors

Solar Cells Lighting

lests

Detection

Instruments

Mechanics Components

Positio

Light Sources

www.toptical.com.tw 02-2346-1510

(standard), dual stage (optional)

ard), P46 (optinal)

dard, 25mm optional

P: 20 x 15mm, 18mm MCP: 14.4 x 10.8mm

nm (quartz, standard), nm (optional)

n the type of the photocathode

n the type of the MCP, up to 35%

e MCP: 4 x 10⁴ MCP: 4 x 10⁶

6 element f/0.8 relay lens on! No vignetting! No pin cushion!



Image intensifier and shutter (schematic)





CCD-Video Unit

CCD Video Chin	Analog Output		Standard	High
	USA, Japan	Elsewhere	Resolution CCD	Resolution CCD
Analog or digital output	analog EIA (RS 170)	analog CCIR	digital 12 or 14 Bit	digital 12 or 14 bit
Resolution (pixel)	768 x 494	752 x 582	780 x 580	1360 x 1024 (1x1) 580 x 512 (2x2) 456 x 342 (ROI)
Pixel size (µm)	8.4 x 9.8	8.6 x 8.3	8.3 x 8.3	4.7 x 4.7
Imaging frequency (analog) Frame rate (digital)	30/60Hz	25/50Hz	12bit/14bit: 33.8/60.8/67.0 fps	12bit/14bit: 10.6/17.9/20.9 fps
Dynamic Range A/D (Bit)	14 Bit, up to 21 Bit (with 4 Spec E spectroscopy software)			
Video Gain	025dB, automatic or manually adjustable through computer RS 232 interface		Full Frame, ROI: 02 Binning: 025dB	20dB
Binning vertical (pixel)	Software		1 (full frame) ,2 (Binr	ning) pixel, ROI
Binning horizontal (pixel)	Software		1 (full frame) ,2 (Binr	ning) pixel, ROI
Image Sensor	ICXAL			
Chip Readout	Correlated double san	npling, dark curre	ent corrected	
Output	$1V_{_{PP}}$ (75 Ω), composit	e video, RS 170/	EIA, CCIR or VGA	
Scan Mode	Field/Frame, selectable through computer RS 232 interface. ICCD camera 4 Picos or 4 Quik E may be Genlocked or supply Fsync-pulse to operate as master clock.			
Gamma	1 or .45, selectable the	rough computer	RS 232 interface	
Internal Synchronization	Free run mode			
External Synchronization	by negative edge TTL pulse (Vinit)			
Cooling of CCD optional	Regulated cooling of of 10 for exposure t	CCD camera uni imes above 100	it to 14°C to minimize Oms. Provides single	dark current by a factor photon sensitivity. No

condensation; eliminates need for vacuum or special nitrogen atmosphere.

4 Quik E control window

Multiple Mode				
Unit D Sequence 0	[u ,50,200,50,200,10	s] _	_	
Warning : Time MCP Gain	exceeds 100 us Voltage 720	P Act	cept 7	
Scan Mode & Frame & Field	Video Acc F Standar CCD Int.	a. d	Frame Gr F None C On Bon	nbber rd
SyncOutput C Vertical C Odd	Trig. Source # Fsync C -Trig	Trigger I C Direct @ STD	Hode Trigger	Gate Fint
Gamma © 0.45 © 1	Video Gain C Auto C Fixed 15	dB	Start Optio	in 17 Au

Shutter control

The advanced, digitally controlled shutter delay feature is the perfect match for your laser, range gating, flow analysis, or many other high speed applications. It is operational in 'scattered light' environments, underwater or for highest speed multi-instrument sequential image acquisition. Multiple direct images with a repetition/ delay time setting as short as 0.3µs can be synchronized with ease to any external TTL source.

internal exposure control

Time (t_s) and delay (t_D) of the gate pulse, or multiple exposure with CPU internally digital programmable	t _s = min t _D =
Trigger propagation delay	<50
Initializing	-Trig
Multiple Exposure	Any

Control of the camera internal Pulse E amplifier via ExtGtP (TTL Pulse) input: Shutter continuous from:	t _s = t, t _D
Trigger propagation delay	<45



www.toptical.com.tw 02-2346-1510 toptical@ms17.hinet.net

	×
	Set Mode
	Single
	Send If
2erl	Connect
	Initialize
	Purt COM1 💌
Þ.	Exit

4 Picos

1.2ns typical* ... 80s,
 a. steps 100ps
 b. ... 80s, min. increments 100ps

Ons, less than 0.01ns jitter

ig, +Trig, or FSync

y sequence, 0.3µs "dead time" between exposures

4 Quik E

= 1.2ns typical* ... DC, $t_D = D_C$ determined by external device

45ns, no jitter

	~
 Imaging 	CD Cameras
Communications	Semiconductors ·
Lighting	Solar Cells ·
Instruments	Tests ·
Sensors	Detection ·
Mechanics · Positioning	Components ·
Light Sources	Lasers ·
2 -	97

Automatic exposure control (optional)

4 Quik E

Pixel by pixel exposure analysis providing automatic light level control by instantaneous 1.2ns typical* ... 15 (18) ms, adjustment of camera shutter speed and Shutter time and MCP gain automatically adjusted, intensifier gain for very wide range of lighting in response to scene illumination. conditions (up to 12 orders of magnitude).

*Actual performance may vary with manufacturer and specimen, please ask for details.



Analog shutter control (schematic)

Mechanical & Environmental Data

Mechanical & Environmental Data, Power Re	equirements
Camera dimensions, without lens (mm, inch)	248 x 110 x 135mm (l x w x h) 9 ³ /4 x 4 ³ /8 x 5 ⁴ /16" (l x w x h)
Camera weight (all in one) (kg / lb)	3kg / 6.6lb
Camera mount (at the bottom plate of the camera)	1/4" x 20 and M8 mounting hole
Operating Humidity (%)	2595%, non condensing
Operating temperature (°C / °F)	0°C – 50°C / 32°F – 122°F
Performance specification	10°C – 40°C / 50°F – 104°F
Operating limits	-10°C – 50°C / 14°F – 122°F
Shock and Vibration	60g accel. shock, 7g Vibration (11 – 200Hz)
Voltage	12 V +5%/-2%





Ż



www.toptical.com.tw 02-2346-1510

Sequence of images for oximetry evaluation © Courtesy M. Crittin et al. IRO, Switzerland 2002



Photgraphy of the retinal oximeter A = ICCD camera *B* = *image* splitter C = fundus camera © Courtesy Courtesy M. Crittin et al. IRO, Switzerland 2002

Multiple 16 bit exposures of Xe flash tube, single discharge, 10 exps during 120µs, continuously swept by scanning mirror.

Spectroscopes	
 Imaging 	CCD Cameras
Communications	Semiconductors ·
Lighting	Solar Cells ·
Instruments	Tests ·
Sensors	Detection ·
Mechanics · Positioning	Components ·
Light Sources	Lasers ·

4 Quik E/... intelligent ICCD Camera Family

digital or analog output			stan	dard		optional
 * Standard: S20 or Super S25 ** Actual performance may vary with manufacturer and specimen, please ask for details. Optional: Bialkali/Fused Silica, Advanced Solar Blind (CsTe), Enhanced S25, Enhanced Blue-UV S20, Wideband S25 WB Two photocathode window materials may be ordered depending upon the required response. 	4 Quik E/01-18**/dig	4 Quik E/01-18**	4 Quik E/01-25**/dig	4 Quik E/01-25**		
Standard Features and Alternatives					, i	
18mm Image Intensifier*, Image area 14.4 x 10.8mm (standard)						
25mm Image Intensifier*, Image area 20 x 15mm						
Minimum Gating Speeds up from 1.2ns typical**						
Integrated, gateable system						
Image Intensifier with Single stage MCP						
Image Intensifier with V-stack dual stage MCP						
Relay Lens Coupling (Image intensifier \rightarrow CCD)						
High Efficiency customized f/0.8 relay lens						
Multiple Exposures						
Shutter dead time 0.3µs						
Integrated Single Trigger Discriminator (STD)						
Analog CCD video output EIA, 768 x 494 pixel or CCIR, 752 x 582 pixel						
SR12, Standard Resolution CCD, 12bit, 780 x 580 pixel						
SR14, Standard Resolution CCD, 14bit, 780 x 580 pixel						
HR12, High resolution CCD, 12bit, 1360 x 1024 pixel						
HR14, High resolution CCD, 14bit, 1360 x 1024 pixel						
Terminal Software and printed manual						
Comfortable case for shipment & storage for free						
Digital Output CameraLink with frame grabber						
Digital output USB 2.0						
Additional Options						
Peltier Cooling						
Special Spectrograph Adapters						
Automatic Exposure Control						
Image Intensifier with Special Photocathodes						
Nikon F-mount Adapter						

Applications

Physical Sciences

Plasma temperature and density analysis Plasma flow analysis Combustion analysis Synchrotron radiation Laser induced fluorescence

Engineering Research

Particle Tracking Velocimetry (PTV) Particle Imaging Velocimetry (PIV) Automative Fuel Injection Spray analysis Wind tunnel studies Stress analysis of ceramics materials

Biological Sciences

Cancer research Fundus imaging spectroscopy X-ray detection Luminescence Time resolved fluorescence

High Speed Imaging

Dynamic Schlieren Phenomena Shock tubes Range gating

Low Light Imaging

Thomson Scattering Raman Spectroscopy Glow Discharge Spectroscopy Semiconductor failure analysis

ICCD - nSec



www.toptical.com.tw 02-2346-1510 toptical@ms17.hinet.net



Specifications are subject to change without notice. No responsibility is assumed for errors or omissions.



Nitrogen-laser-based system for oral cancer diagonis developed at Center for Advanced Technology, Indore, India © with courtesy of CAT, Indore, India

Spectroscopes	
 Imaging 	CCD Cameras
Communications	Semiconductors ·
Lighting	Solar Cells •
Instruments	Tests ·
Sensors	Detection ·
Mechanics · Positioning	Components ·
Light Sources	Lasers ·
۰ ۱	404



XX Rapid Frame

Fastest multi-framing ICCD camera system



Features

- Three or Four intensified CCD channels (4 Picos or 4 Quik E)
- Total of 3/6 or 4/8 frames from one trigger
- Fastest shutter down to 0.2ns (based on 4 Picos ICCD cameras) in single mode
- Jitter <10ps
- Ultra fast recording of up to 4 full frame resolution images with 0.01ns interframing time
- Double shutter/frame and multiple exposures (each module)
- · Perfectly spectral flat image splitter
- Single Photon detection
- Spectral Sensitivity of photocathodes from UV NIR (220 1300nm)
- USB2.0 or analog output
- High Dynamic Range: 14 Bit (>36 Bit with 4 Spec E spectroscopy software)
- · Effortless Image/Data storage and retrieval
- Internal or external trigger
- Free Terminal software "ready-to-use"

SPECIFICATIONS

Unique Features	XXRF based on 4 Picos	XXRF based on 4 Quik E
Shortest gating time each module	0.2ns	1.2 ns typically*
Jitter	< 0.01ns	< 0.01ns
Multiple exposures, "dead" time between exposures	any sequence 0.3µs	
Gate repetition rate	3.3MHz burst, 200kHz continue	ous (each)
Intensifier output coupling	customized distortion free f/0.8 Image Splitter	relay lens
Image Splitter	Spectrally flat mirror system	

* Minimum gate time 1.2ns typically (depends on manufacturer and specimen)

Special emphasis was laid on the utmost precision of the timing control of the system

The XXRapidFrame ultra high speed digital imaging system opens a new range of timing capabilities. The timing precision (iitter) is four times more accurate than that of dedicated delay generators (e.g. DG535) the minimum gate/delay times and gate/delay time steps are 100 times shorter than previous "state of the art" systems.

Exact reproducible digital setting of delay and exposure times is standard.

Very high system integration permits compact physical size of the XXRapidFrame camera. Distortion free imaging due to advanced proximity focused MCP (Micro Channel Plate) image intensifier, a optimized custom designed f0.8 relay coupling lens, and the use of highest quality CCD array for best sensitivity and resolution results in superior total optical performance.

14 Bit High Dynamic Range can be increased with 4 Spec E Spectroscopy PC Software up to >36 Bit (11 orders of magnitude) with multiple frames integrated.



Digital or Analog Data

Image Splitting Concepts

The wide spread semitransparent image splitting setup does not provide equal spectral distribution to both channels. Our spectrally flat Mirror Image Splitter does have identical amount of light of all wavelengths to all channels.

Semi transparent image splitter



The "conventional" method used to split an image into two identical secondary images of both lower intensity uses a cubic or pellicle beam splitter. It leads to a non predictable intensity ratio because the original image is usually not monochromatic. Therefore the photon flux at every single pixel is determined by the ratio of the integrals of the spectral distribution weighted with the R/T ratio of the beam splitter. Any wavelength change in between frames, e.g. as a function of temperature, will modify the splitting ratio. Such resulting images cannot be numerically corrected





CCD Cameras

Semiconductors munications

Solar Cells

lests

Detection

Compo

Instruments

Mechanics

Positio

_ight Sources

Imaging

Sol





Mirror image splitter



Using totally reflective mirrors for all sub images yields the intensity division in between the subs being only a function of the geometry of the mirrors. This approach is easily extended into the UV spectral region. A filter placed in each channel's optical path to produce e.g. a three-color image of an experiment. With different delay times in the various channels it could also be used to retrieve 3-D spatial information.

Image Intensifier

Image Intensifier specifications

Image intensifier type (proximity focused MCP)	single stage (standard), dual stage (optional)
Phosphor material	P43, P46
Image intensifier diameter (mm)	18mm, 25mm
Image area of the relay lens	18mm MCP: 14.4 x 10.8mm, 25mm MCP: 20 x 15mm
Wavelength range, subject to window design	165 – 920nm (quartz, standard), 165 – 1300nm (optional)
Spectral Sensitivity of MCP (nm)	2201300nm, depends on the type of the photocathode
Quantum Efficiency (Q.E.) (see curves below)	depends on the type of the MCP, up to 35%
Gain (4k steps) (V_{MCP} =1000V) control via RS 232 digital setup	single stage MCP: 4×10^4 dual stage MCP: 4×10^6
Signal to noise (db @ µLx)	46dB min @ 0.5µLx
Coupling phosphor (MCP \rightarrow CCD)	customized 6 element f/0.8 relay lens No distortion! No vignetting! No pin cushion!

Spectral Sensitivity of Photocathodes (Wavelength in nm)

Standard 18 mm				
S20UV	В	approx. 165 - 820nm		
S25 IR (Super S25)	н	approx. 350 - 920nm		
Ор	otiona	l 18 mm		
S20 UV (MgF2)	A	approx. 110 - 820nm		
Broadband	J	approx. 190 - 920nm		
			3	
			Enha	
			N	

Deviations of up to ±25% from the above typical spectral sensitivity curves are possible. The position of the curves can vary ±20nm. The input window material limits the spectral response of the photocathode in the shorter wavelengths. The window materials and their transmission limits are: quartz (165nm), MgF2 (110nm).



Image intensifier and shutter (schematic)





www.toptical.com.tw 02-2346-1510





CCD-Video Unit

	Analog	J Output	Progressive	Standard	High
CCD Video Chip	USA, Japan	Elsewhere	Scan CCD	Resolution CCD	Resolution CCD
Analog or digital output	analog EIA	analog CCIR	analog, VGA 30/60Hz or 60/120Hz	digital 12 or 14 Bit	digital 12 or 14 bit
Resolution (pixel)	768 x 494	752 x 582	640 x 480	780 x 580	1 3 6 0 x 1 0 2 4 (1x1) 580 x 512 (2x2) 456 x 342 (ROI)
Pixel size (µm)	8.4 x 9.8	8.6 x 8.3	9.8 x 9.8	8.6 x 8.6	4.8 x 4.8
lmaging frequency (analog) Frame rate (digital)	30/60Hz	25/50Hz	30/60 (30/60Hz) 60/120/200/240/350 (60/120Hz)	12bit/14bit: 33.8/ 60.8/ 67.0 fps	12bit/14bit: 10.6/17.9/20.9 fps
Video Gain	025dB, autor computer RS 2	matic or manuall 232 interface	y adjustable through	Full Frame, R Binning: 025	OI: 020dB dB
Binning vertical (pixel)	Software			1 (full frame) , 2 (Binning) pix	kel, ROI
Binning horizontal (pixel)	Software			1 (full frame) , 2 (Binning) pix	kel, ROI
Output	1V _{PP} (75 Ω), c VGA	omposite video,	RS 170/EIA, CCIR or	USB 2.0	
Dynamic Range A/D (Bit)	14 Bit, up to 27	1 Bit (with 4 Spec	c E spectroscopy softwa	are)	
Chip Readout	Correlated dou	uble sampling, da	ark current corrected		
Internal Synchronization	Free run mode				
External Synchronization	by negative edge TTL pulse (Vinit)				
Signal to noise	46dB min @ 0.5µLx				
Cooling of CCD optional Regulated cooling of CCD camera unit to 14°C to minimize dark current by a factor of for exposure times above 100ms. Provides single photon sensitivity. No condensation				ent by a factor of 10 . No condensation;	

for exposure times above 100ms. Provides single photon sensitivity. No condensation; eliminates need for vacuum or special nitrogen atmosphere.

Mechanical & Environmental Data, Power Requirements

XXRapidFrame (based on 4 Picos or 4 Quik E)	Three
Camera dimensions, without lens (I x w x h)	480 x 300
Camera weight (all in one) (kg / lb)	25kg / 55l
Camera mount (at the bottom plate of the camera)	3/8" x 20 a
Operating Humidity (%)	2595%,
Operating temperature (°C / °F)	0°C – 50°
Performance specification	10°C – 40
Operating limits	-10°C – 5
Shock and Vibration	60g accel
Voltage	12 V +5%

Shutter control

The advanced, digitally controlled shutter delay feature is the perfect match for your laser, range gating, flow analysis, or many other high speed applications. It is operational in 'scattered light' environments, underwater or for highest speed multi-instrument sequential image acquisition. Multiple direct images with a repetition/ delay time setting as short as 0.3µs can be synchronized with ease to any external TTL source.

Internal exposure control	XXRapi
Fime (t_s) and delay (t_D) of the gate pulse, or multiple exposure with CPU internally ligital programmable	t_s = 0.2ns min. steps 1 t_p = 0 80s increments
rigger propagation delay	<65ns, less
nitializing	-Trig, +Trig,
Aultiple Exposure	Any sequen

* mimimum gate time 1.2ns typically, depends on manufacturer and specimen.

Ź



· Imaging CCD Cameras

> Communications Semiconductors

Solar Cells Lighting

Tests · Instruments

Detection Sensors

Mechanics Components

Positioning

www.toptical.com.tw 02-2346-1510

e Channels

Four Channels

x 350 mm 650 x 350 x 380mm

30kg / 66lb

and M8 mounting hole

non condensing

°C / 32°F – 122°F

0°C / 50°F – 104°F

50°C / 14°F – 122°F

I. shock, 7g Vibration (11 – 200Hz)

6/-2%

idFrame P

XXRapidFrame E

80s, 10ps , min. 10ps

 $t_s = 1.2ns^*$ typically... 80s, min. steps 100ps $t_{\rm D}$ = 0 ... 80s, min. increments 0.1ns

than 10ps jitter

or FSync

ice, 0.3µs "dead time" between exposures

Lasers Light Sources

External exposure contro	I
--------------------------	---

XXRapidFrame P

device

XXRapidFrame E

 $t_s = 1.2ns^*$ typically. ∞ , t_D , ∞ t, t_D determined by external device

Control of the camera internal Pulse E amplifier via ExtGtP (TTL Pulse) input: Shutter continuous from:

Trigger propagation delay

<45ns, no jitter

t_s = 0.2ns ... ∞, t_D, ∞

t, t_{D} determined by external



Analog shutter control (schematic)

Applications

ICCD - Rapid Frame

Physical Sciences

Plasma temperature and density analysis Plasma flow analysis Combustion analysis Synchrotron radiation Laser induced fluorescence

Biological Sciences

Cancer research Fundus imaging spectroscopy X-ray detection Luminescence Time resolved fluorescence



High Speed Imaging

Dynamic Schlieren Phenomena Shock tubes Range gating

Low Light Imaging

Thomson Scattering Raman Spectroscopy Glow Discharge Spectroscopy Semiconductor failure analysis

XXRapid Frame

digital or analog output

			stand	ard	optional	
*25mm Standard: S20 or S25, minimum gating time 1.5 ns typically (depends on manufacturer and specimen) **18mm Standard: S20 or Super S25, minimum gating time 1.2 ns typically (depends on manufacturer and specimen) Optional: Bialkali/Fused Silica, Advanced Solar Blind (CsTe), Enhanced S25, Enhanced Blue-UV S20, Wideband S25 WB Two photocathode window materials may be ordered depending upon the required response	XXRF-E (4QuikE)	XXRF-Edig (4QuikE/dig	XXRF-P (4Picos)	XXRF-Pdig (4Picos/dig)		Spectroscopes
Standard						maging
Cating Speeds from 1 2ps /1 5ps typically					<u> </u>	
Gating Speeds from 0.2ns $\rightarrow \infty$	-	-				
Gating Speeds from 0.215 $\rightarrow \infty$						Comn
						nunic
						ation
Image Intensifier with Vistack dual stage MCP						0
						Light
						ing
Multiple Exposures						
Shutter dead time 0.3us						
Integrated Single Trigger Discriminator (STD)						nstru
Analog CCD video output EIA, 768 x 494 pixel or CCIR, 752 x 582 pixel	•	_		-		uments
Progressive Scan CCD, VGA, 640 x 480 pixel						
Standard Resolution CCD, 12bit, 750 x 580 pixel						Sens
Standard Resolution CCD, 14bit, 750 x 580 pixel						sors
High resolution CCD, 12bit, 1360 x 1024 pixel						
High resolution CCD, 14bit, 1360 x 1024 pixel						
USB 2.0						Mech
Terminal Software and printed manual						anics
Comfortable aluminum case for storage for free						·
Additional Options			· · · ·		···	osition
Peltier Cooling						ing
Special Spectrograph Adapters						
Image Intensifier with Special Photocathodes						Ligh
Nikon F-mount Adapter						it Sou

t



Quantum Leap

The modular high speed image intensifier module continuous or gated operation (e.g. with Time Gain Module, optional)



Features

- · Compact module with built-in image intensifier optimized relay optics
- · Gated or continuous operation
- Gated System: Gating times from 0.2 DC (with TGE) or 1.2ns or (with TGN)
- 18 or 25 mm Image intensifier
- Spectral Sensitivity of photocathode from UV NIR (110 1300nm)
- · Very compact high voltage power supply and pulse amplifier
- Photonic gain up to 106, 0 1000V, 10-turn potentiometer or remotely controlled
- Single Photon detection
- Improves image contrast and S/N ratio
- · Customized 6-element distortion free lens coupling between image intensifier and CCD
- Distortion free imaging
- Needs only 12V, 300mA electrical power
- · Precision adjustment rings for back focusing to video camera

SPECIFICATIONS

Unique Features	Quantum Leap	Quantum Leap N Shortest gating	/E (gated)
Shortest gating time	non gateable	1.2 ns typically*	E: 0.2ns
Remote control	not applicable	RS 232 (standard)	
Optical input	c-mount (standard), N	likon F-mount (optional)
Optical input window	quartz		
Sensitivity corresponding to conventional film	$1 \times 10^{12} \text{fc} = 50 \times 5^7 \text{AS}$	SA	
Multiple exposures, "dead" time between exposures	any sequence 0.3µm		
Gate repetition rate	3.3MHz burst, 200kHz	z continuous	
Intensifier output coupling	customized distortion	free f/0.8 relay lens	
CE certified	yes		

Programmable control parameter entry via RS 232 (with Time & Gain Module), or 10-turn potentiometer Distortion free imaging due to advanced proximity focused MCP (Micro Channel Plate) image intensifier and use of highest quality CCD array for best sensitivity and resolution.



Image Intensifier

Image Intensifier specifications Image intensifier type (proximity focused MCP) Phosphor material Image intensifier diameter (mm) Image area of the relay lens Wavelength range, subject to window design Spectral Sensitivity of MCP (nm) Quantum Efficiency (Q.E.) (see curves below) Gain (4k steps) (V_{MCP} =..1000V) control via RS 232 digital setup

Coupling phosphor $(MCP \rightarrow CCD)$



Spectral Sensitivity of Photocathodes (Wavelength in nm)						
Standard 25 mm		Optional 25 mm				
S20	С	approx. 165 - 820nm	S20 UV(MgF2)	А	approx. 110 - 820nm	
S25	F	approx. 200 - 840nm	Solar Blind (CsTe)	G	approx. 180 - 340nm	
Standard 18 mm		Bialkali		approx. 165 - 600nm		
S20UV	В	approx. 165 - 820nm	Enhanced S25 (glass)	D	approx. 165 - 820nm	
S25 IR (Super S25)	н	approx. 350 - 920nm	Enhanced S25 (glass)	I	approx. 270 - 900nm	
Optional 18 mm		Wideband S25 WB	к	approx. 200 - 900nm		
Broadband	J	approx. 190 - 920nm	S1	Е	approx. 700 - 1300nm	

Deviations of up to ±25% from the above typical spectral sensitivity curves are possible. The position of the curves can vary ±20nm. The input window material limits the spectral response of the photocathode in the shorter wavelengths. The window materials and their transmission limits are: quartz (165nm), MgF2 (110nm).





match M

Shutter control (with optional Time & Gain Module)

The advanced, digitally controlled shutter delay feature is the perfect match for your laser, range gating, flow analysis, or many other high speed applications. It is operational in 'scattered light' environments, underwater or for highest speed multi-instrument sequential image acquisition. Multiple direct images with a repetition/ delay time setting as short as 0.3µs can be synchronized with ease to any external TTL source.

Internal exposure control	Qu
Time (t_s) and delay (t_D) of the gate pulse, or multiple exposure with CPU internally digital programmable	$t_s = 0.2ns$ min. step $t_D = 0$ increme
Trigger propagation delay	<65ns, le
Initializing	-Trig, +T
Multiple Exposure	Any sequ

External exposure control	Quantum Le
Control of the camera internal Pulse E amplifier via ExtGtP (TTL Pulse) input: Shutter continuous from:	$t_s = 0.2ns \dots \infty, t_s$ t, t_p determined device
Trigger propagation delay	<45ns, no jitter





Mechanical & Environmental Data

Applications

Physical Sciences

Plasma flow analysis Combustion analysis

Synchrotron radiation

Quantum Leap, modular screw on Module

Mechanical & Environmenta	I Data, Power Requirements	
mensions (I x b x h), without lens	60 x 236 x 150mm, 2 3/8 x 5 3/8 x 5 7/8"	
veight (all in one) (kg / lb)	1.7kg	*25mm Standard: S20 or S25
mount (at the bottom plate of the camera)	1/4" x 20 and M8 mounting hole	**18mm Standard: S20 or Super S25 Optional: Bialkali/Fused Silica, Advanced Solar Blind (CsTe),
ting Humidity (%)	2595%, non condensing	Enhanced S25, Enhanced Blue-UV S20, Wideband S25 WB
ting temperature (°C / °F)	0°C – 50°C / 32°F – 122°F	Two photocathode window materials may be ordered depending upon the required response.
rmance specification	10°C – 40°C / 50°F – 104°F	
ating limits	-10°C – 50°C / 14°F – 122°F	
k and Vibration	60g accel. shock, 7g Vibration (11 – 200Hz)	Standard
ge	12 V +5%/-2%	Gated system
		Gating Speeds from 1.2ns $\rightarrow \infty$ (18mm MCP)
		Gating Speeds from 0.2ns $\rightarrow \infty$
ICATIONS		Image Intensifier with 18mm MCP
ical Sciences		Image Intensifier with 25mm MCP
a tomporature and density analysis		Image Intensifier with Single stage MCP
a flow analysis		Image Intensifier with V-stack dual stage MCP
ustion analysis		Lens Coupling
induced fluorescence		High Efficiency customized relay lens
		Two c-Mount outputs
neering Research		Programmable control parameter via RS232
le Tracking Velocimetry (PTV)		Multiple Exposures
le Imaging Velocimetry (PIV)	High Speed Imaging	Shutter dead time 0.3µs
ative Fuel Injection analysis	Dynamic Schlieren Phenomena	Terminal Software and printed manual
unnel studies	Shock tubes	Comfortable case for shipment & storage for free
analysis of ceramics materials	Range gating	Additional Options
aical Sciences	Low Light Imaging	Nikon F-mount Adapter (for 25mm MCP)
		Image Intensifier with Special Photocathodes

Specifications are subject to change without notice. No responsibility is assumed for errors or omissions.

Engineering Research

Laser induced fluorescence

Particle Tracking Velocimetry (PTV) Particle Imaging Velocimetry (PIV) Automative Fuel Injection Spray analysis Wind tunnel studies Stress analysis of ceramics materials

Biological Sciences

Fundus imaging spectroscopy X-ray detection Luminescence Time resolved fluorescence

Thomson Scattering Raman Spectroscopy Glow Discharge Spectroscopy Semiconductor failure analysis







4 Spec E

Spectra Acquisition Software for pulsed and time resolved Spectroscopy



4 Spec E is a software application for the PCs with Pentium or subsequent microprocessor. It is designed for use with a frame grabber and any video camera for video imaging spectroscopy, where image data is grouped into bands that run either parallel or perpendicular to the video scan lines. Intensity data from each band may be plotted as intensity graphs, as pseudo-color or grayscale 16- bit/pixel images, or as a superimposition of image and intensity plots.

4 Spec E is designed to run under Windows 9x, ME or NT/2000/XP or subsequent on the IBMPC. Analysis of stored curve data highresolution plots as well as grayscale or full-color images can be printed to any IBM-compatible printer. Graphics can alternatively be achieved to a TIFF file for export to drawing or publishing programs. Data may also be exported as text or binary integers for entry into spreadsheets, image processing software, etc. The Data can also be output as real-time video into the Image Pro Sequence files (*.SEQ).

4 Spec E also supports Dynamic Data Exchange (DDE) protocol as a DDE server for curve data.

Intensity data from an image may be grouped and summed for plotting into a maximum of 8000 different 16 bit per pixel linear arrays, called curve memories. Curve length is adjustable up to 8000 pixels. 4 Spec E lets you choose the specific image rows to include in each curve, and supports arithmetic processing of curves and images.

4 Spec E offers a full range of camera-control and data-acquisition functions. Data for the full frame grabber image or a grouped image array can be acquired and viewed as grayscale or pseudo-color images, or as curves in a variety of 2D and 3D plot modes. Any type of image or plot can be repetitively scanned in live mode, yielding convenient, real-time data presentation. Real-time feature include automatic darkbackground image subtraction, automatic image flat-fielding, and live display of image pixel position and intensity.

4 Spec E is designed for real-time data acquisition with a PC or compatible computer incorporating a Pentium, or subsequent processor. 4 Spec E communicates with the installed hardware via a MS-Windows software driver which is included on the same disk. These drivers are available for a growing number of frame grabbers and slow scan data acquisition systems. If you wish to use a device which is not yet supported, please let us know. Not only are we continually adding additional drivers, we can also provide special drivers upon

Specifications

Sensitivity	more than 1 count/photoelectron/ pixel. Up to 80 seconds integration time on CCD detector
Image configuration	1 line scan to 582 lines horizontal and vertical binning
Dynamic Range (with cooled detector)	More than 21 Bit (2 Mio:1),all lines integrated (binned), dynamic expansion active.
Display options	Stores up to 30,000 curve memories in RAM, 32 Bit/pixel curve memories; 256 level gray scale or 256 level pseudocolor; 2-D and 3-D overlays; peak finder; wave length, wave number, and photon energy calibration image function for laser profile displays.

4Spec - Samples.crv (10 eurves)

Bie Edit Set Ejecute Style Window Process Setup (mgl Cgl Help

Import (Image Mode)	Edit (Curve Mode)	24
aw Dete to Instage	Qeda Dir+U	General Serlings .
entFile to Conve(t) - entFile to CAL Darve -	Deter Clear gl Deters All Ab-Det	Arto Deckground Subtock Arto Bell Peld Conect
aller Fieldette ,	Smooth Curice Citing	Auto Trendore to Absorbe Auto Trendore to Logarithm
Nort (Curve/Catalog Mode)	Converte Picer Offic	Satas Background Cave Satas Rehiground Cave
HERD CALLS:	Pasar Plot Layrad	Satas Ato Braik Curve
	Edit (Catalog Mode)	Evenine.
per. Col-0 eve Col-S	Out Shar-Out Dray Oni-tes	Groete Dverlay Plot Deuley Catalog
eon i	Essle Shifting Dirice Del Celeni An-Del	Pain Scel Soquerce Acquire One Exposure Acquire One Exposure Brocksport
da Serv póes hi Fie	Edit (Image Mode)	PagastGuastrage
nça Sahal	Dear	Acquire Rost-time Sumon -
nt Al Davan	ClearAl	Disaley Active Image Disaley Destigravital Image Disaley Flet Flet I made in age
Export (Image Mode)		
Invega es Peur Date Invega es TET File . Invega es DAP File		International States of St
oort (Curve/Catalog Mode)		FatFieldCorectinope
Curve(s) as Text		Fig Inego Hatonia Fig Inego Vetoni
Property Rev Determination of the second sec		



· Imaging

nications

Solar Cells

_ 🗆 ×

www.toptical.com.tw 02-2346-1510

Data Acquisition

From an 8-bit per pixel frame grabber image, rows and columns are selected for reading inti a 16 bit/ pixel computer image having 65.000:1 dynamic range. The exposure sequence is programmed by choosing the number of image frames (from one to many thousands) to be accumulated into 32bit/ pixel curve memories for plotting. Any exposure may be externally triggered.

The user can include automatic static- or dynamic background subtraction, flat-field correction, and conversion to logarithmic or absorbance scales.

Depending upon computer speed and image size, up to 15 - 30 frames per second may be captured.

Acquire Real-time series . This command enables the scanning of a series of frames at video rate. All images, images segments, curves, etc. may be automatically saved to disk or RAM disk in real time.

	Style (Curve Mode)	244	= -
	Like Color	Curve Memories .	N.
DHY	FBColer .	Sites Burway.	
orkoe	Control of the second s	interallianticat Chini	B
Sec.	Zura V Okcin Od+7	Scel Searche Oti+G	<u>e</u>
dence	Calibrate Areas	Curve Definition OX+N	lts
	The Onl	Memory We de	
	+ Auss	Add Sara	
OHF	Contraction of the second second	Seve Calena	
- Cal-6		Lated Options	S S S
	Style (Image Mode)		SU
	C. Der Cres	thelp	0
		Contents	S
#3	Image Color +	About	
Chi-M	0.48.1999034	Ļ	
CrinB	and the second s	SAME NED	
Chief		stat porer per later	
nt CrixD		4 SPEC Help Index	\leq
C19-2	Carve Enods or Image	Menu Commanda	ec
Cost		the Annu Contracts	ha
n. FI		Last Manual Control of State	nic.
CrisA	Positive Ott-1	Lineals, Hera, Carttania,	<i></i> ы
	Negative Chi+2 «	Mitchen Mena-Laterards	•
	Pseudocolor (NH)	THE MOULDMENTS	Po
		Control Palaties	Sit .
	Same Const	Fairlie in Carlos Mode	0
-	Representation Chief	Franklin in annalise Marine	nin
See.	Calibration	How to	Q
wan Cei-B	Arregister -	over the surgicity of the West Name	
e	Live Worker . Ote-V	Linked, Assumption in Caller Linke Free Freed Contractions	
	Camera Control	inter Sectors of Links	
	Established Galage		<u> </u>
		Appendix	Int
		Training	S
		The second contracts	0

Versatile and Adaptable

Plotting

Two- and three dimensional overlay plots of curves may be generated with a single mouse click. The 3D plots have handles for dragging to resize the image and change the perspective. In addition, a Channel Profile plot which graphs intensity along a slice parallel to the Z axis may be added. Various plot choices available include line, bar scatter, and cityscape styles, with optional color fill.

Displays

Images may be shown in real time with 256 levels of gray or in 256 level pseudo color and may be viewed in actual size or resized to fit the display window. The option, contour Image, transforms curve data into patterned, pseudo color image. In addition, 4 Spec E features a live laser-beam profiler display.

Focusing

Adjustments to focus are simply done by using live-image or repeat-scan curve plot displays. The user is able to monitor the intensity of any pixel in real time in any image or plot with the moveable cursor. 4 Spec E employs a peak finder feature which can automatically monitor the positions, widths and areas of peaks in a curve during scanning.

Super Zoom

Super Zoom uses the mouse to zoom and pan the curve plot displays in camera live action (even in 3D). This allows the user to zoom in and out and pan the display while scanning in real-time with the camera.

Curve Catalog

This feature makes use of Curve Icons to give a single window view of all curve memories. A curve icon may be plotted by double-clicking a curve icon. To erase drag to 'curve trash' and click. Curve icons may be cut, copied, and pasted to edit memories.

	🚪 4Spec - SAMPLE	S.CRV (10 curves)			_ O X
	Ello Edit Set Ex	ecute Style <u>Windo</u>	w Process Setup	Ord Impl Help	
	1: Sample_1	2: Sample_2	3: Sample_3	4: Sample_4	5: Sample_5
	سالعانات أساس	بالغالل والم	ملياللساس	باللاسي	
, Э	6: Sample_6	7: Sample_7	8: Sample_B	9: Sample_9	10: Sample_10
	•	Scale		Croats Ovorley	J

Automatic Spike Filter

The automatic spike filter removes high energy events (electrical discharges and other interferences) fro images acquired during long exposure.

Export curve data

Export curve data may be in ASCII text format. Images may be exported in 16-bit binary or TIFF formats. Images may be imported from TIFF, or binary-integer-data files.

Calibration

Spectra may be calibrated in wavelength, Raman cm-1 shift or other units with linear or cubic-splinefits. Peak positions and widths are located with the 'intelligent' peak finder. Spectral sensitivity may be calibrated with a tungsten lamp.



Acquire Real-time series

This command enables the scanning of a series of frames at video rate. The frames are stored in the frame grabber buffer. Selected frames can then be transferred into the active image.

	Execute	Sty <u>l</u> e	<u>W</u> indow	<u>P</u> rocess	<u>S</u>
	Run Me				
	Create	F3			
	Display	Ctrl+M			
	Run So	Ctrl+R			
	Acquire	e One E	xposure	Ctrl+E	
	Acquire	Ctrl+D			
	Repea	Ctrl+Z			
	Repea	Ctrl+L			
÷	Acquire	F2			
	Display	y Active	Image	Ctrl+A	
	Display	y Backg	ground Imag	ge	
	Display	y Flat-Fi	ield Image		

ICCD





CCD Cameras

Imaging

8

nications

Lighting

Instruments

Sel Detection

Mechar Components

lics

Positio

ight Sources

lotors

Solar Cells

lests

www.toptical.com.tw 02-2346-1510



Movies and live scans

may be run as animated 3D waterfall graphs which dramatically show spatial or time evolution of spectral data.

Plots

Plots may be printed on any PC compatible printer, including color ink-jet printers. Plots may be pasted to the clipboard for copying into other documents.

Frame grabber

Drivers for the following frame grabbers are available: MV-Gamma (4 Quik E/4 Picos, analog), MV-Titan (XXRapid Frame) and .he MV-CL CameraLink for digital ICCD cameras are available.

USB 2.0

For digital ICCD cameras with USB 2.0 no frame grabber is needed.